

ElectroForce™ Apex 1



Site Preparation Guide

Table of Contents

| | |
|----------------------------------|-------|
| Table of Contents | 2 |
| Ideal Setup | 3 |
| System Components | 4 |
| Instrument Measurements | 5-7 |
| Instrument..... | 5 |
| Power Supply | 6 |
| Controller | 7 |
| Utility Requirements | 8-9 |
| Power | 8 |
| Gas | 9 |
| Accessories | 10-15 |
| Hot/Cold Chamber..... | 10-11 |
| Environmental Bath | 12-13 |
| ACS-3..... | 14-15 |
| Site Preparation Checklist | 16 |
| TA Instrument Offices | 17 |



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab




Customer

Ideal Setup



IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor and ceiling space and a rigid laboratory bench that is level and is in a vibration-free environment. Bench must be rated to support several hundred pounds.



Distance from the wall:
15.24 cm (0.5 ft) min.

Allow additional 86.36 cm
(34 in) from the wall for
optional hot/cold chamber

Table width: 152.4 cm (5 ft)


Table depth: 121.92 cm (4 ft)

Floor space: 91.44 cm x
91.44 cm (3 ft x 3 ft)

System Components



MAIN SYSTEM COMPONENTS



- | | |
|--------------------|------------------|
| A Test Instrument | D Power Supply |
| B Computer Tower | E Controller |
| C Computer Monitor | F Emergency Stop |

Instrument Measurements



APEX 1 – FLOOR STANDING FRAME



Max Height: 247.65 cm (97.5 in)
Min Height: 203.2 cm (80 in)

Width: 68.58 cm (27 in)

Depth: 68.58 cm (27.0 in)

Weight (Axial only): 257 kg (567 lbs)
Weight (Axial-Torsion): 290 kg (639 lbs)



APEX 1 – BENCHTOP

Max Height: 179.1 cm (70.5 in)
Min Height: 134.6 cm (53 in)

Width: 62 cm (27 in)

Depth: 53.3 cm (21 in)

Weight: 201 kg (443 lbs)



Instrument Measurements




AXIAL POWER SUPPLY FOR APEX 1

Height: 40.64 cm (16 in)


Width: 38.1 cm (15 in)

Depth: 40.64 cm (16 in)

Weight: 21.8 kg (48 lbs)



TORSION POWER SUPPLY FOR APEX 1



Height: 30.48 cm (12 in)

Width: 33 cm (13 in)

Depth: 40.64 cm (16 in)

Weight: 14.3 kg (31.5 lbs)

Instrument Measurements




CONTROLLER FOR APEX 1

Height: 11.43 cm (4.5 in)

Width: 38 cm (15 in)

Depth: 40.64 cm (16 in)


Weight: 4.1 kg (9 lbs)



Utility Requirements



POWER

| Item | Requirement |
|----------------------|---|
| Power Supply | <ul style="list-style-type: none">• Axial: 208–230V, 50–60 Hz, 10A• Torsion: 208–230V, 50–60 Hz, 2A• Neutral to Ground (NG) voltage max 0.5 volt• Safety ground per local regulation• Single phase |
| Controller Power | <ul style="list-style-type: none">• 207–230V, 50–60 Hz Hz, 1.2A |
| Power cords provided | <ul style="list-style-type: none">• 6-20P plug for 230V systems in North America and Japan• International: Line power cord provided is based on country  6-20P |



Use power cords with plugs appropriate for your circuit.



Supply voltages lower than indicated may result in a degradation of performance.



Ensure that the mains assigned do not also supply power to noise generating equipment nearby, such as motors, welders, transformers, etc.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab




Customer


Utility Requirements



GAS

| Item | Requirement | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|----------------------|--|----------------------|-------------------|-------------------|-------------------|------------------|--------------------------------------|--|-----|--|--|--|-------|--------|--------------------------------------|--|------------------|----------------|--------------|-------------------|-------------------|------------------|--------------------------|---|----------|-------|------|---|-----------------|---|--------------|---|-----------|---------|-------|---|-----------------|---|------------|-----------------|---|----------|---------|---|----------------|---|----------|-------------------|---|---|----------|---|---------------|---|---------|---|---|---|-----------|---|---------------|---|---|---|---|---|---|-----|----------------|---|---|
| System gas | Air | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Pressure | 552–621 kPa (80–90 psig) NOTE: Air gauge provided with system | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Flow | 45.3 L/min (1.6 CFM) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Connections | ¼-inch push-to-connect or ¼-inch NPT female to filter/dryer/regulator with ¼-inch push-to-connect removed ¹ | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Conditions | <table><tr><th rowspan="2">ISO8573-1:2010 Class</th><th colspan="3">Solid Particulate</th><th rowspan="2">Concentration</th><th colspan="2">Water</th><th>Oil</th></tr><tr><th colspan="3">Maximum number of particles per m³</th><th>Vapor</th><th>Liquid</th><th>Total oil (aerosol, liquid, & vapor)</th></tr><tr><th></th><th>0.1 – 0.5 micron</th><th>0.5 – 1 micron</th><th>1 – 5 micron</th><th>mg/m³</th><th>Pressure dewpoint</th><th>g/m³</th><th>ppm (mg/m³)</th></tr><tr><td>1</td><td>≤ 20,000</td><td>≤ 400</td><td>≤ 10</td><td>-</td><td>≤ -94°F (-70°C)</td><td>-</td><td>0.008 (0.01)</td></tr><tr><td>2</td><td>≤ 400,000</td><td>≤ 6,000</td><td>≤ 100</td><td>-</td><td>≤ -40°F (-40°C)</td><td>-</td><td>0.08 (0.1)</td></tr><tr><td>Recommended → 3</td><td>-</td><td>≤ 90,000</td><td>≤ 1,000</td><td>-</td><td>≤ -4°F (-20°C)</td><td>-</td><td>0.83 (1)</td></tr><tr><td>Minimum Specs → 4</td><td>-</td><td>-</td><td>≤ 10,000</td><td>-</td><td>≤ +37°F (3°C)</td><td>-</td><td>4.2 (5)</td></tr><tr><td>5</td><td>-</td><td>-</td><td>≤ 100,000</td><td>-</td><td>≤ +45°F (7°C)</td><td>-</td><td>-</td></tr><tr><td>6</td><td>-</td><td>-</td><td>-</td><td>≤ 5</td><td>≤ +50°F (10°C)</td><td>-</td><td>-</td></tr></table> | ISO8573-1:2010 Class | Solid Particulate | | | Concentration | Water | | Oil | Maximum number of particles per m ³ | | | Vapor | Liquid | Total oil (aerosol, liquid, & vapor) | | 0.1 – 0.5 micron | 0.5 – 1 micron | 1 – 5 micron | mg/m ³ | Pressure dewpoint | g/m ³ | ppm (mg/m ³) | 1 | ≤ 20,000 | ≤ 400 | ≤ 10 | - | ≤ -94°F (-70°C) | - | 0.008 (0.01) | 2 | ≤ 400,000 | ≤ 6,000 | ≤ 100 | - | ≤ -40°F (-40°C) | - | 0.08 (0.1) | Recommended → 3 | - | ≤ 90,000 | ≤ 1,000 | - | ≤ -4°F (-20°C) | - | 0.83 (1) | Minimum Specs → 4 | - | - | ≤ 10,000 | - | ≤ +37°F (3°C) | - | 4.2 (5) | 5 | - | - | ≤ 100,000 | - | ≤ +45°F (7°C) | - | - | 6 | - | - | - | ≤ 5 | ≤ +50°F (10°C) | - | - |
| ISO8573-1:2010 Class | Solid Particulate | | | Concentration | Water | | Oil | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | Maximum number of particles per m ³ | | | | Vapor | Liquid | Total oil (aerosol, liquid, & vapor) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 0.1 – 0.5 micron | 0.5 – 1 micron | 1 – 5 micron | mg/m ³ | Pressure dewpoint | g/m ³ | ppm (mg/m ³) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | ≤ 20,000 | ≤ 400 | ≤ 10 | - | ≤ -94°F (-70°C) | - | 0.008 (0.01) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | ≤ 400,000 | ≤ 6,000 | ≤ 100 | - | ≤ -40°F (-40°C) | - | 0.08 (0.1) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Recommended → 3 | - | ≤ 90,000 | ≤ 1,000 | - | ≤ -4°F (-20°C) | - | 0.83 (1) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Minimum Specs → 4 | - | - | ≤ 10,000 | - | ≤ +37°F (3°C) | - | 4.2 (5) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | - | - | ≤ 100,000 | - | ≤ +45°F (7°C) | - | - | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | - | - | - | ≤ 5 | ≤ +50°F (10°C) | - | - | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |





¹ Interconnection between the system and the supplied filter/regulator equipment is via ¼-inch push-to-connect fittings; use supplied tubing.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab




Customer

Accessories



HOT/COLD CHAMBER MEASUREMENTS



Height: 30.48 cm (12 in)

Width: 33 cm (13 in)

Depth: 40.64 cm (16 in)

Weight: 14.3 kg (31.5 lbs)



HOT/COLD CHAMBER CONTROLLER MEASUREMENTS

Height: 14 cm (5.5 in)

Width: 45.2 cm (17 in)

Depth: 35 cm (13 in)

Weight: 2.26 kg (5 lbs)



Accessories



HOT/COLD CHAMBER REQUIREMENTS

| Item | Requirement |
|-------------|--------------------------------------|
| Gas | Liquid nitrogen |
| Pressure | 152–345 kPa (22–50 psig) |
| Connections | ½ -inch SAE, 45 degree flare fitting |
| Power | 230V, 50–60 Hz Hz,9.6A/2200W, 1 ph |



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab




Customer

Accessories



ENVIRONMENTAL BATH MEASUREMENTS



Height: 30.48 cm (12 in)

Width: 33 cm (13 in)

Depth: 40.64 cm (16 in)

Weight: 14.3 kg (31.5 lbs)



ENVIRONMENTAL BATH CONTROLLER MEASUREMENTS

Height: 14 cm (5.5 in)

Width: 36.58 cm (14 in)

Depth: 27 cm (10 in)

Weight: 1.36 kg (3 lbs)



Accessories



ENVIRONMENTAL BATH REQUIREMENTS

| Item | Requirement |
|-------|-----------------------------|
| Power | 207–230 VAC, 50/60 Hz, 3.0A |



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Accessories



AIR CHILLER SYSTEM (ACS-3) MEASUREMENTS



Height: 112 cm (44 in)

Width: 37 cm (14.5 in) WITHOUT Chiller Panel

Width: 56 cm (22 in) WITH Chiller Panel

Depth: 56 cm (22 in)

Weight: 121 kg (267 lbs) WITHOUT Chiller Panel

Weight: 128 kg (282 lbs) WITH Chiller Panel




CHILLER PANEL MEASUREMENTS

Height: 64 cm (25 in)

Width: 56 cm (22 in)

Depth: 38.1 cm (15 in)

Weight: 7.25 (16 lbs)



Accessories



AIR CHILLER SYSTEM (ACS-3) REQUIREMENTS

Requirements



- 50 Hz: 207–252 VAC (refer to the serial number plate on the rear of the unit. The ACS is line frequency specific.)
- 60 Hz: 216–252 VAC (refer to the serial number plate on the rear of the unit. The ACS is line frequency specific.)
- 8A, 9A, 9.1A, or 11A (refer to the serial number plate on the rear of the unit)
- US sites require an L6-20 single-phase outlet



NEMA L6-20 plug



- Gas:
 - Air or nitrogen
 - Pressure: 6.9 bar (100 psig)
 - Flow rate: 200 SLPM
 - Temperature: 20–30°C
 - Dew point: -40°C (-40°F) **NOTE:** Dew point is specified at operating pressure. Supplying dryer air at a lower dew point will extend continuous operation.



- Lab Environment (must be below 25°C): 21°C–24°C = Ideal
- Leave 20 cm (8 in) of space in the front and back of the ACS for ventilation



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Site Preparation Checklist



ElectroForce Apex 1



Enough bench/floor space for instrument and computer

- ☐ Table width: 1.5 m (5 ft)
- ☐ Table depth: 1.2 m (4 ft)
- ☐ Floor space (for floor standing frame): 91.44 cm x 91.44 cm (3 ft x 3 ft)
- ☐ If optional Hot/Cold Chamber purchased, allow an additional 87.8 cm (34 in) from the wall so the system can slide back when not in use



Power supply is:

- ☐ Axial: 208–230V, 50–60 Hz, 10A
- ☐ Torsion: 208–230V, 50–60 Hz, 2A

Controller power is:

- ☐ 207–230V, 50–60 Hz, 1.2A



Purge gas

- ☐ Air: 552–621 kPa (80–90 psig)
- ☐ Flow rate: 45.3 L/min (1.6 CFM)
- ☐ 5 μ m filtered
- ☐ Dew point \leq 10°C



Optional Hot/Cold Chamber:

- ☐ Gas: Liquid nitrogen
- ☐ Pressure: 152–345 kPa (22–50 psig)



Optional ACS-3:

- ☐ Gas: Air or nitrogen
- ☐ Pressure: 6.9 bar (100 psig)
- ☐ Flow rate: 200 SLMP
- ☐ Dew point: -40°C (-40°F)
- ☐ L6-20 single phase outlet for US sites



- ☐ The Customer assumes responsibility for any damage that occurs when the instrument is moved by someone other than a trained TA Instruments Service Representative.

I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.

If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.

Customer DD / MM / YYYY

Company City State Country

Please send a signed copy of the completed checklist to your local Service representative.

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com